SWISS FFP2

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SWISS FFP2 NR RESPIRATOR
N°ART. 1602 E | EN 149 : 2001 + A1 : 2009

WWW.SWISS-FFP2.CH

The foldable SWISS FFP2 respirator offers high wearing comfort and effective protection of the respiratory tract against solid and liquid aerosols according to EN 149: 2001 + A1: 2009.

INSTRUCTIONS FOR USE

Before using the mask, the following instructions must be read and followed.

Failure to follow the instructions for use can result in illness or death!

- The user of the mask must be instructed on how to use it correctly.
- Masks do not protect against gases and vapors.
- The oxygen content in the breathing air must be minimum 19.5% by volume.
- Respiratory protection must not be used if the concentration, type or characteristics of the harmful substances are not known.
- The mask must be changed immediately if it is damaged or the breathing resistance increases.
- To ensure the safety of the wearer, the mask must be changed regularly

- (at least every 4 hours) and little manipulated. It is not reusable.
- Never make changes or modifications to the mask.
- Check the condition of the mask before use.
- The mask is not flame resistant and is not intended for use in potentially explosive atmospheres.

INSTRUCTIONS FOR PUTTING ON THE MASK

Wash your hands before touching the mask.

Unfold the mask, especially the nose clip area and pull on the elastic straps, making sure that the foldable band is facing upwards.

Form the bridge of the nose before positioning the mask on the face. Insert the chin first, then place the mask on the nose. Then attach the elastic bands.

Shape the clip over the nose area to the contour of the face.

Check the correct positioning of the mask: Place both hands on the mask and exhale vigorously. No air should escape at the edges.

Warning: Facial hair such as beards can affect the required firm fit of the mask and can cause leaks.

STORAGE AND SHELF-LIFE

Shelf life when properly stored max. 3 years.

Check expiry date on masks before use. Storage and transport clean and dry between -5°C to + 25°C at max. 80% humidity in original packaging. Protect against UV radiation!

DISPOSAL

In the event of damage, high pollution or if breathing resistance is significantly increased, the mask should be changed and disposed of in accordance with regional regulations.

CONFORMITY ASSESSMENT

The particle filtering half mask meets the requirements of the European Directive (EU) 2016/425. Type examination and conformity assessment procedures were carried out with the participation of the following notified body:

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Download the declaration of conformity from: www.pro-geneve.ch/SwissFFP2/Declaration_de_conformite

INFO

We are happy to help you choose the right respirator. The selection of the most suitable personal protective equipment (PPE) depends on the respective situation and may only be carried out by competent persons. For questions and further information please contact PRO Foundation. www.pro-geneve.ch

We thank you for your trust, with the SWISS FFP2 mask you have chosen a Swiss quality product, made by PRO. PRO is a private non-profit foundation recognized as being of the public utility. PRO's mission is to professionally reintegrate by offering a job in a suitable environment, people excluded from the traditional economy, due, for example, to a disability. PRO staff are trained in the production of FFP2 masks and constant quality control. PRO manufactures high quality respiratory protection masks.

Products are 100% made in Switzerland and meet the highest quality standards according to CE guidelines. The company has its own internal test laboratory, whereby the outstanding filter quality is constantly controlled and ensured.

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